

**Chris Stavroudis solution published in WAAC Newsletter, January 2013:**

**(Thanks to Elise Rousseau and her post to the CIPP list, David Goist, Mary-Lou Florian, Hilary Kaplan, Jane Bassett, Rustin Levinson and Marc Williams for their thoughts and contributions.)**

**For surfaces that can be exposed to a small amount of bleach, eg., the reverse of paintings , stretchers, frames, etc.**

**47 ml 100% isopropanol  
23 ml 100% ethanol  
7 ml 3% hydrogen peroxide  
23 ml distilled water**

**The same recipe can be made from materials available from any well stocked pharmacy as follows:**

**44 ml 91% isopropanol rubbing alcohol  
30 ml 70% ethyl alcohol rubbing alcohol  
7 ml 3% hydrogen peroxide  
19 ml distilled water**